



**K100-0104.5**  
**4.5 - 10 pH Test Strips**  
**Instructions for Use**

**Indications for Use**

K100-0104.5, 4.5 - 10 pH Test Strips are indicated for use in testing the pH of water used for preparing dialysate. There are no other indications for use.

**Directions**

- 1) Remove one strip from the container and immediately recap the container.
- 2) Dip the strip into the water to be tested and move back and forth until the color on the indicator pads has stopped changing.
- 3) Remove the strip and shake off excess liquid. Compare the strip to the color chart on the container.
- 4) Match the strip as closely as possible to one of the color scales on the container. Read the pH listed above the matched colors. Colors that are between two scales indicate a pH that is between the two values for those scales.

**Warnings and Cautions**

- Do not touch the indicator pads.
- Keep all unused test strips in the original container.
- Replace cap immediately/tightly after removing a strip.
- Do not use test strips from an opened or unopened container after expiration date printed on the container label.
- Do not allow the test strip to come in contact with liquids or with work surfaces that may contaminate the strip.
- Do not leave test strips in areas exposed to vapors of any type.
- *Caution:* When used as a medical device, Federal Law restricts this device to sale by or on order of a physician.

**Interferences**

No interferences are known.

**Storage**

Keep container cool and dry.

Store test strips at a temperature between 39° F (4°C) and 86° F (30°C).