



K100-0104.6.3
6.0 – 7.7 pH Test Strips
Instructions for Use

Indications for Use

K100-0104.6.3, 6.0 – 7.7 pH Test Strips are indicated for use in testing the pH of water used for preparing dialysate. There are no other indications for use.

Directions

- 1) Remove one strip from the container and immediately recap the container.
- 2) Dip the strip into the water to be tested and move back and forth until the color on the indicator pad has stopped changing.
- 3) Remove the strip and shake off excess liquid. Compare the strip to the color chart on the container.
- 4) Match the strip as closely as possible to one of the color blocks on the container. Read the pH listed above the matched color. A color that is between two scales indicates a pH that is between the two values for those scales.

Warnings and Cautions

- Do not touch the indicator pad.
- Keep all unused test strips in the original container.
- Replace cap immediately/tightly after removing a strip.
- Do not use test strips from an opened or unopened container after expiration date printed on the container label.
- Do not allow the test strip to come in contact with liquids or with work surfaces that may contaminate the strip.
- Do not leave test strips in areas exposed to vapors of any type.
- *Caution:* When used as a medical device, Federal Law restricts this device to sale by or on order of a physician.

Interferences

No interferences are known.

Storage

Keep container cool and dry.

Store test strips at a temperature between 39° F (4°C) and 86° F (30°C).