



K100-0117 E-Z Chek® 6.8 - 8.5 pH Test Strips Instructions for Use

Indications for Use

K100-0117 E-Z Chek® 6.8 - 8.5 pH Test Strips are indicated for use in testing the pH of acid/bicarbonate dialysate and bicarbonate concentrate. There are no other indications for use.

Directions for Acid/Bicarbonate Dialysate or Bicarbonate Concentrate

- 1) Remove one strip from the bottle and immediately recap the bottle.
- 2) Dip the strip into the solution to be tested and move back and forth for one second.
- 3) Remove the strip and shake off excess liquid. After 10 seconds, compare the strip to the color chart on the bottle.
- 4) Match the strip as closely as possible to one of the color scales on the bottle. Read the pH listed above the matched colors. Colors that are between two scales indicate a pH that is between the two values for those scales.

Warnings and Cautions

- Do not touch the indicator pads.
- Keep all unused test strips in the original bottle.
- Replace cap immediately/tightly after removing a strip.
- Do not use test strips from an opened or unopened bottle after expiration date printed on the bottle label.
- Do not allow the test strip to come in contact with liquids or with work surfaces that may contaminate the strip.
- Do not leave test strips in areas exposed to vapors of any type.
- *Caution:* When used as a medical device, Federal Law restricts this device to sale by or on order of a physician.

Interferences

No interferences are known.

Storage

Keep container cool and dry.

Store test strips at a temperature between 39° F (4°C) and 86° F (30°C).