

Urinalysis

Combi 10[®] SGL
Optimized urine test strip

**CLIA-
waived**



Best available vitamin C protection

+ Safe, early diagnosis of serious diseases

Desiccant enclosed in the cap

+ Convenient and safe

Premium color charts ideally match reaction colors


+ Clear read-off, reliable results

Distributed by:




(800) 647-3873

• Medi-Test Combi 10[®] SGL



Best available vitamin C protection

Comparative data clearly show that urine test strips Medi-Test Combi 10[®] SGL have the best available protection of the blood and glucose test against influences of vitamin C. With other test strips, vitamin C can prevent the oxidation reaction used to detect blood and glucose in urine to cause false negative readings. Using Medi-Test Combi 10[®] SGL **blood and glucose** in urine **will be safely detected** and patients can receive treatment at the earliest stage with optimal chances for success.



Desiccant enclosed in cap

Test strips Medi-Test Combi 10[®] SGL come in tubes with a desiccant containing cap. The desiccant ensures a low humidity in the tube and the longevity of the test strips. Many other packages contain a simple desiccant pack that may easily be spoiled or accidentally removed. The desiccant in the cap is **convenient and safe**. Closing the cap safely protects the test strips.



Suitable for visual and instrumental evaluation

The reaction colors on the test strip exactly match the colors on the label. This is ensured by a high quality printing process using specially mixed full colors. The evaluation of Medi-Test Combi 10[®] SGL is always **clear and reliable** providing perfectly accurate results. Medi-Test Combi 10[®] SGL urine test strips can also be evaluated using the CLIA-waived strip reader URYXXON[®] Relax. Instrumental readings eliminate the subjectivity of visual color interpretation and minimize risks associated with manual transcription.



Ordering information

Item#	Description
K100-0125	Combi 10 [®] SGL Reagent test strips for blood, urobilinogen, bilirubin, protein, nitrite, ketones, glucose, pH, specific gravity and leucocytes in urine by the dip and read technique. Tube of 100 strips
K100-0130	URYXXON [®] Relax for the instrumental evaluation of urine test strips Combi 10 [®] SGL.

Distributed by:



(800) 647-3873



KATEN100078
FlyerCombi10SGL USA1/3/0/01.2008 PD
Printed in Germany

Medi-Test Combi 10[®] SGL