



K100-0104 pH Test Strips - Instructions for Use

- 1) Remove one strip from the container and immediately recap the container.**
- 2) Immerse the strip into the solution to be tested and move the strip back and forth in the solution.**
- 3) If the color change on the strip color pads has stopped changing after 5 seconds, then compare the strip to the scales on the container. Match the strip as closely as possible to one of the container scales and read the pH listed above the matched colors. Colors that are between two scales indicate a pH that is between the two values for those scales.**

Note: Allow the test strip to stay immersed in the solution, until the color no longer changes, before making a comparison to the scales on the container (for some test solutions this may be longer than 5 seconds). The pH-Fix test strips will be suitable for use up to the expiration date stated on the box, even when they are outside the box, as long as they are stored dry, at room temperature, and protected from light.

RPC • Tucson, AZ • Minneapolis, MN • USA • (520) 888-5551 • Fax (520) 888-5557

Product Quality at Discount PricingSM



K100-0104.5 pH Test Strips - Instructions for Use

- 1) Remove one strip from the container and immediately recap the container.**
- 2) Immerse the strip into the solution to be tested and move the strip back and forth in the solution.**
- 3) If the color change on the strip color pads has stopped changing after 5 seconds, then compare the strip to the scales on the container. Match the strip as closely as possible to one of the container scales and read the pH listed above the matched colors. Colors that are between two scales indicate a pH that is between the two values for those scales.**

Note: Allow the test strip to stay immersed in the solution, until the color no longer changes, before making a comparison to the scales on the container (for some test solutions this may be longer than 5 seconds). The pH-Fix test strips will be suitable for use up to the expiration date stated on the box, even when they are outside the box, as long as they are stored dry, at room temperature, and protected from light.

RPC • Tucson, AZ • Minneapolis, MN • USA • (520) 888-5551 • Fax (520) 888-5557

Product Quality at Discount PricingSM