



Tube Occluding Clamp/ Lightweight/Low Cost Hemostat



Key Features

- Plastic clamp costs only a fraction of its metal counterpart
- Used to occlude medical grade tubing from 1 mm to 8 mm.
- Smooth jaws won't damage tubing
- Jaws remain aligned when tubing is occluded
- Reusable, easy to use, lightweight, safe, and effective
- Available in Dark Blue, Sky Blue, White, Green, Red, and Yellow

Materials

- Composed of lightweight, durable nylon (nylon 6,6)
- Aluminum rivet (Won't rust when cleaning)

Cleaning

- Clamps can be cleansed using warm water and a detergent, mixed in a sonic cleaner. Diluted bleach solutions are acceptable, however it can cause discoloration of the plastic. Cold disinfectant solutions, such as Cidex, can be used to disinfect the clamps during a soak. Whenever possible, follow medical device cleaning protocols in your facility as well as warning labels on cleaning solutions to prevent or minimize damage during the process. Do not use strong acid-based cleaners.